



eDNA Monitoring Research

- USGS round-robin assay evaluation
- USGS streamgauge eDNA sampling

Adam Sepulveda (asepulveda@usgs.gov)

USGS Northern Rocky Mountain Science Center

USGS Upper Midwest Environmental Science Center


USGS Forest and Rangeland Ecosystem Science Center

USGS Western Fisheries Research Center



Need

1. How to sample?
2. Where to sample?



Round-robin evaluation of eDNA assays

eDNA is powerful early detection tool but...

Uncertainty associated with eDNA results

- Indirect method that detects DNA
 - DNA presence \neq Species presence
- Often impossible to corroborate with direct methods
- Newer science with checkered early history



Towards reducing uncertainty of eDNA results...

Identify lab-based protocols that provide repeatable & reproducible results.

- Lab environment can be totally controlled
- Repeatable: same sample & conditions = same results
- Reproducible: same sample, different conditions same result
- Cornerstone good science & diagnosis
 - confidence that initial results not obtained by chance or errors.



Approach, 5 assays & 4 labs

- **Public, probe-based quantitative PCR assays**
 1. DRE16S: zebra & quagga
 2. DRE2: zebra
 3. ZebCOI: zebra
 4. ZebCYT: zebra
 5. DREQM: quagga

- **Four USGS labs**
 - Boise, ID; Bozeman, MT; LaCrosse, WI; Seattle, WA



Approach

With mussels

- Seneca Lake (ZM & QM)
- Lake Michigan (ZM & QM)
- San Justo Res (ZM)
- Lake Mead (QM)




No mussels, spiked with ZM + QM DNA

- Yellowstone River
- Jackson Lake
- Columbia River



Approach

- Partners collected samples using SOP
- Double-blind
 - Sample prep
 - Sample analysis
- Identical samples
- 2832 samples, 8497 qPCR reps



Repeatability ~ 90%

Likelihood that two identical samples analyzed under the same conditions (same assay or same lab) will yield same result

- 1 = results are the same

Assay	Mean	CI95
DRE16S	0.92	0.03
DREQM	0.93	0.02
DRE2	0.85	0.03
ZEBCOI	0.91	0.03
ZEBCYT	0.93	0.03
Lab	Mean	CI95
1	0.93	0.03
2	0.92	0.03
3	0.92	0.03
4	0.86	0.04



Reproducibility ~ 0.90

Likelihood that two identical samples analyzed under different conditions (assay or lab) will yield same result

- 1 = results are the same

	Mean	CI95	
Assays	0.92	0.01	
Labs	0.92	0.01	
Assay	DRE16S	ZEBCOI	ZEBCYT
DRE2	0.90	0.90	0.90
ZEBCYT	0.93	0.93	
ZEBCOI	0.93		
DREQM	0.95		
Lab	4	3	2
1	0.89	0.93	0.94
2	0.90	0.93	
3	0.87		




Conclusions & Next Steps

- High probability (~90%) that identical samples analyzed under the same or different assays or labs will yield the same result
- Dreissenid mussel eDNA: increase confidence / decrease uncertainty in results
- Benchmark for new assays and /or labs
- Next steps: evaluate eDNA field sampling protocols




USGS Streamgauge eDNA sampling

- 10,000+ USGS streamgages
- USGS hydro techs visit ~ monthly
- Incorporate eDNA sampling into visits
 - Cheap and effective way to help *Safeguard the West?*



USGS streamgages below reservoirs

- Reservoirs
 - invasion hotspots
 - critical infrastructure
- USGS gages often located below reservoirs
- eDNA at gages
 - Mussel (-), vulnerable
 - Mussel (+)



Objectives

- 1. eDNA surveillance for dreissenid mussels**
- 2. Additional cost of incorporating eDNA into USGS streamgauge program**
- 3. Are streamgages adequately located for bio-surveillance of upstream reservoirs?**
 - What influences increase detection rates?



Approach

- **Hydro-techs trained in SOP**
 - 3 – 10 water samples per visit
 - May – September 2018
 - Shipped/analyzed at USGS LaCrosse, WI
 - **Target species**
 - **Dreissenid mussels**
 - Likely all no-detects... hard to learn
 - **Kokanee salmon & Yellow perch**
 - Reservoir 'obligates' at varying abundance
 - Used these results to learn
-



Results

1. Dreissenid mussel DNA

- Mussel pos. waters: all samples positive
- Mussel neg. waters: no sample met criteria for being scored as positive

2. Additional costs:

- \$400 – 500 for three samples at a streamgage
- Additional 2 hr of hydro-tech time


3. Streamgages adequately located

- Kokanee and yellow perch DNA: yes
-




Results: detection rate in 1 sample

Kokanee



Yellow perch




Results: how many samples?

Kokanee

- diminishing returns > 6

Yellow perch

- Less certain, but likely > 10



Take homes

Streamgage eDNA sampling

- always detected mussel DNA when mussels are present
- sensitive even when fish at low abundance
- cost-effective

Next steps

- What are we missing by limiting to streamgage?
 - Send USGS hydro techs to reservoirs rather than, or in addition to, streamgage?



"Preliminary Information-Subject to Revision. Not for Citation or Distribution."

